

Technical Bulletin

Guidelines for Measuring Lead in Blood: Repeat Analysis and Confirmation Testing¹



Magellan
Diagnostics

Guideline 13: Repeat Analysis of the Original Specimen

1. If initial result is $\geq 5 \mu\text{g/dL}$, reanalyze, if sample volume permits, to rule out contamination.
 - If there is insufficient volume, e.g. capillary specimen, report initial result and refer patient for confirmatory testing. (See Guideline 14)
 - If a large discrepancy is encountered, the specimen should be analyzed a 3rd time, discarding the outlier and reporting the average of the remaining two.***
2. For any result exceeding $5 \mu\text{g/dL}$, refer the patient for confirmatory testing.

Guideline 14: Confirmatory Testing

1. When blood lead results $\geq 5 \mu\text{g/dL}$ are obtained, the laboratory must either:
 - Refer the specimen to a CLIA certified lab for confirmatory testing or
 - If there is insufficient sample for another analysis, see Guideline 15.
 - Preliminary results may be released with a comment that 'results of the confirmatory test are pending'.

Guideline 15: Reporting POC Results to Clinical Health Care Providers

1. Blood lead results $\geq 5 \mu\text{g/dL}$ require further investigation.
 - Consider using a comment similar to: Blood lead levels $\geq 5 \mu\text{g/dL}$ indicate that a child has been exposed to lead at levels higher than most children. The blood lead level should be confirmed using a venous blood sample. Since no safe blood lead level has been identified, no detectable level should be considered 'normal'.

***Large Discrepancies are defined as:

Blood Lead Result	Discrepancy Limit
5 - 20 $\mu\text{g/dL}$	> 3 $\mu\text{g/dL}$
21 - 40 $\mu\text{g/dL}$	> 4 $\mu\text{g/dL}$
> 40 $\mu\text{g/dL}$	> 10%

Reference:

1. *Guidelines for Measuring Lead in Blood Using Point of Care Instruments, Draft Report of the Laboratory Workgroup of the Advisory Committee on Childhood Lead Poisoning Prevention Of the Centers for Disease Control and Prevention October 24, 2013*
http://www.cdc.gov/nceh/lead/ACCLPP/20131024_POCguidelines_final.pdf [accessed 24 Jan 2014].